

ADAMTS12 Antibody



PACO24763

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Metalloprotease that may play a role in the degradation of COMP. Cleaves also alpha-2 macroglobulin and aggrecan. Has anti-tumorigenic properties.

Gene ID:

ADAMTS12

Uniprot

P58397

Synonyms:

A disintegrin and metalloproteinase with thrombospondin motifs 12 (ADAM-TS 12) (ADAM-TS12) (ADAMTS-12) (EC 3.4.24. -), ADAMTS12


Immunogen:

Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 12 protein (26-229AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO24763 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO24763 at dilution of 1:100.