

ADAMTS4 Antibody, Biotin conjugated



PACO24770

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves aggrecan, a cartilage proteoglycan, and may be involved in its turnover. May play an important role in the destruction of aggrecan in arthritic diseases. Could also be a critical factor in the exacerbation of neurodegeneration in Alzheimer disease. Cleaves aggrecan at the '392-Glu-|-Ala-393' site.

Gene ID:

ADAMTS4

Uniprot

O75173

Synonyms:

A disintegrin and metalloproteinase with thrombospondin motifs 4 (ADAM-TS 4) (ADAM-TS4) (ADAMTS-4) (EC 3.4.24.82) (ADMP-1) (Aggrecanase-1), ADAMTS4, KIAA0688

Immunogen:

Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 4 protein (213-319AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A