

ADORA1 Antibody, Biotin conjugated



PACO24794

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor for adenosine. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Gene ID:

ADORA1

Uniprot

P30542

Synonyms:

Adenosine receptor A1, ADORA1

Immunogen:

Recombinant Human Adenosine receptor A1 protein (202-235AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A