
Product Information

Size:

50ug

Reactivity:

Mouse, Rat, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Binds estrogens, fatty acid, and metals.

Gene ID:

Afp

Uniprot

P02772

Synonyms:

Alpha-fetoprotein (Alpha-1-fetoprotein) (Alpha-fetoglobulin), Afp


Immunogen:

Recombinant Mouse Alpha-fetoprotein protein (18-605AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Afp antibody at 10 μ g/ml

Lane 1: HepG2 whole cell lysate

Lane 2: Rat liver tissue

Lane 3: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 68 kDa

Observed band size: 68 kDa

.