

Aif1 Antibody, FITC conjugated



PACO24828

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Actin-binding protein that enhances membrane ruffling and RAC activation. Enhances the actin-bundling activity of LCP1. Binds calcium. Plays a role in RAC signaling and in phagocytosis. May play an role in macrophage activation and function. Promotes the proliferation of vascular smooth muscle cells and of T-lymphocytes. Enhances lymphocyte migration. Plays a role in vascular inflammation.

Gene ID:

Aif1

Uniprot

P55009

Synonyms:

Allograft inflammatory factor 1 (AIF-1) (Ionized calcium-binding adapter molecule 1) (MRF-1) (Microglia response factor), Aif1, Iba1 Mrf1

Immunogen:

Recombinant Rat Allograft inflammatory factor 1 protein (2-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A