

ALAS2 Antibody, Biotin conjugated



PACO24854

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, mitochondrial matrix, mitochondrion, 5-aminolevulinate synthase activity, coenzyme binding, glycine binding, cellular iron ion homeostasis, erythrocyte differentiation, heme biosynthetic process, hemoglobin biosynthetic process.

Gene ID:

ALAS2

Uniprot

P22557

Synonyms:

5-aminolevulinate synthase, erythroid-specific, mitochondrial (ALAS-E) (EC 2.3.1.37) (5-aminolevulinic acid, synthase 2) (Delta-ALA synthase 2) (Delta-aminolevulinate synthase 2), ALAS2, ALASE ASB

Immunogen:

Recombinant Human 5-aminolevulinate synthase, erythroid-specific, mitochondrial protein (50-587AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A