

ALDH5A1 Antibody, HRP conjugated



PACO24877

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes one step in the degradation of the inhibitory neurotransmitter gamma-aminobutyric acid, (GABA).

Gene ID:

ALDH5A1

Uniprot

P51649

Synonyms:

Succinate-semialdehyde dehydrogenase, mitochondrial (EC 1.2.1.24) (Aldehyde dehydrogenase family 5 member A1) (NAD(+)-dependent succinic semialdehyde dehydrogenase), ALDH5A1, SSADH

Immunogen:

Recombinant Human Succinate-semialdehyde dehydrogenase, mitochondrial protein (48-535AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A