

Annexin A5 Antibody



PACO24932

Product Information

Size:

50ug

Reactivity:

Cynops pyrrhogaster

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Gene ID:**Uniprot**

P70075

Synonyms:

Annexin A5 (Annexin V) (Annexin-5)

Immunogen:

Recombinant Cynops pyrrhogaster Annexin A5 protein (1-323AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A