

# APAF1 Antibody, HRP conjugated



PACO24965

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Oligomeric Apaf-1 mediates the cytochrome c-dependent autocatalytic activation of pro-caspase-9 (Apaf-3), leading to the activation of caspase-3 and apoptosis. This activation requires ATP. Isoform 6 is less effective in inducing apoptosis.

**Gene ID:**

APAF1

**Uniprot**

O14727

**Synonyms:**

Apoptotic protease-activating factor 1 (APAF-1), APAF1, KIAA0413

**Immunogen:**

Recombinant Human Apoptotic protease-activating factor 1 protein (104-415AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A