

PACO25032

Product Information

Size:

50ug

Reactivity:

Guinea Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acid, , hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

Gene ID:

ALB

Uniprot

H0V003

Synonyms:

Serum albumin (allergen Can f 3), ALB

Immunogen:

Serum Albumin purified from Guinea Pig serum (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A