

# ARHGAP25 Antibody, HRP conjugated



PACO25037

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

**Gene ID:**

ARHGAP25

**Uniprot**

P42331

**Synonyms:**

Rho GTPase-activating protein 25 (Rho-type GTPase-activating protein 25), ARHGAP25, KIAA0053

**Immunogen:**

Recombinant Human Rho GTPase-activating protein 25 protein (17-280AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A