

ARHGAP25 Antibody, Biotin conjugated



PACO25039

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GTPase activator for the Rho-type GTPases by converting them to an inactive GDP-bound state.

Gene ID:

ARHGAP25

Uniprot

P42331

Synonyms:

Rho GTPase-activating protein 25 (Rho-type GTPase-activating protein 25), ARHGAP25, KIAA0053

Immunogen:

Recombinant Human Rho GTPase-activating protein 25 protein (17-280AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A