

CENPB Antibody, FITC conjugated



PACO25082

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Interacts with centromeric heterochromatin in chromosomes and binds to a specific 17 bp subset of alphoid satellite DNA, called the CENP-B box. May organize arrays of centromere satellite DNA into a higher-order structure which then directs centromere formation and kinetochore assembly in mammalian chromosomes.

Gene ID:

CENPB

Uniprot

P07199

Synonyms:

Major centromere autoantigen B (Centromere protein B) (CENP-B), CENPB

Immunogen:

Recombinant Human Major centromere autoantigen B protein (1-165AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A