

# CENPB Antibody, Biotin conjugated



PACO25083

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Interacts with centromeric heterochromatin in chromosomes and binds to a specific 17 bp subset of alphoid satellite DNA, called the CENP-B box. May organize arrays of centromere satellite DNA into a higher-order structure which then directs centromere formation and kinetochore assembly in mammalian chromosomes.

**Gene ID:**

CENPB

**Uniprot**

P07199

**Synonyms:**

Major centromere autoantigen B (Centromere protein B) (CENP-B), CENPB

**Immunogen:**

Recombinant Human Major centromere autoantigen B protein (1-165AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A