

# Atoh1 Antibody, HRP conjugated



PACO25121

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional regulator. Activates E box-dependent transcription in collaboration with TCF3/E47, but the activity is completely antagonized by the negative regulator of neurogenesis HES1. Plays a role in the differentiation of subsets of neural cells by activating E box-dependent transcription.

**Gene ID:**

Atoh1

**Uniprot**

P48985

**Synonyms:**

Protein atonal homolog 1 (Helix-loop-helix protein mATH-1) (mATH1), Atoh1, Ath1

**Immunogen:**

Recombinant Mouse Protein atonal homolog 1 protein (1-351AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A