

ATP12A Antibody, Biotin conjugated



PACO25127

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the hydrolysis of ATP coupled with the exchange of H(+) and K(+) ions across the plasma membrane. Responsible for potassium absorption in various tissues.

Gene ID:

ATP12A

Uniprot

P54707

Synonyms:

Potassium-transporting ATPase alpha chain 2 (EC 3.6.3.10) (Non-gastric H(+)/K(+) ATPase subunit alpha) (Proton pump), ATP12A, ATP1AL1

Immunogen:

Recombinant Human Potassium-transporting ATPase α chain 2 protein (171-306AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A