

AZIN1 Antibody, Biotin conjugated



PACO25179

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Antizyme inhibitor protein that positively regulates ornithine decarboxylase (ODC) activity and polyamine uptake by counteracting the negative effect of antizyme OAZ1, OAZ2 and OAZ3 on ODC1 activity. Inhibits antizyme-dependent ODC degradation by binding to antizymes. Releases ODC1 from its inactive complex with antizymes, leading to formation of the catalytically active ODC1.

Gene ID:

AZIN1

Uniprot

O14977

Synonyms:

Antizyme inhibitor 1 (AZI) (AZI1) (Ornithine decarboxylase antizyme inhibitor), AZIN1, OAZI OAZIN

Immunogen:

Recombinant Human Antizyme inhibitor 1 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A