

BCAT1 Antibody, FITC conjugated



PACO25218

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the first reaction in the catabolism of the essential branched chain amino acid, leucine, isoleucine, and valine.

Gene ID:

BCAT1

Uniprot

P54687

Synonyms:

Branched-chain-amino-acid, aminotransferase, cytosolic (BCAT(c)) (EC 2.6.1.42) (Protein ECA39), BCAT1, BCT1 ECA39

Immunogen:

Recombinant Human Branched-chain-amino-acid, aminotransferase, cytosolic protein (4-276AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A