

BCAT2 Antibody, HRP conjugated



PACO25221

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the first reaction in the catabolism of the essential branched chain amino acid, leucine, isoleucine, and valine. May also function as a transporter of branched chain alpha-keto acid, .

Gene ID:

BCAT2

Uniprot

O15382

Synonyms:

Branched-chain-amino-acid, aminotransferase, mitochondrial (BCAT(m)) (EC 2.6.1.42) (Placental protein 18) (PP18), BCAT2, BCATM BCT2 ECA40

Immunogen:

Recombinant Human Branched-chain-amino-acid, aminotransferase, mitochondrial protein (28-392AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A