

PACO25232

---

## Product Information

**Size:**

50ug

**Reactivity:**

Canine

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acid, , hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.

**Gene ID:**

ALB

**Uniprot**

P49822

**Synonyms:**

Serum albumin (allergen Can f 3), ALB

**Immunogen:**

Serum Albumin purified from Canine serum (Native Protein).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A