

# Bcl2 Antibody, FITC conjugated



PACO25238

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Suppresses apoptosis in a variety of cell systems including factor-dependent lymphohematopoietic and neural cells. Regulates cell death by controlling the mitochondrial membrane permeability. Appears to function in a feedback loop system with caspases. Inhibits caspase activity either by preventing the release of cytochrome c from the mitochondria and/or by binding to the apoptosis-activating factor (APAF-1).

**Gene ID:**

Bcl2

**Uniprot**

P10417

**Synonyms:**

Apoptosis regulator Bcl-2, Bcl2, Bcl-2

**Immunogen:**

Recombinant Mouse Apoptosis regulator Bcl-2 protein (5-205AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A