

# BDH1 Antibody, Biotin conjugated



PACO25247

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

D-beta-Hydroxybutyrate dehydrogenase mitochondria has been purified to apparent homogeneity. The membrane-bound enzyme is first released by phospholipase A digestion of the mitochondria.

**Gene ID:**

BDH1

**Uniprot**

Q02338

**Synonyms:**

D-beta-hydroxybutyrate dehydrogenase, mitochondrial (EC 1.1.1.30) (3-hydroxybutyrate dehydrogenase) (BDH) (Short chain dehydrogenase/reductase family 9C member 1), BDH1, BDH SDR9C1

**Immunogen:**

Recombinant Human D-beta-hydroxybutyrate dehydrogenase, mitochondrial protein (47-343AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A