

Bglap Antibody, FITC conjugated



PACO25254

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Constitutes 1-2% of the total bone protein. It binds strongly to apatite and calcium.

Gene ID:

Bglap

Uniprot

P04640

Synonyms:

Osteocalcin (Bone Gla protein) (BGP) (Gamma-carboxyglutamic acid, containing protein), Bglap, Bglap2

Immunogen:

Recombinant Rat Osteocalcin protein (50-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A