

LGB Antibody



PACO25284

Product Information

Size:

50ug

Reactivity:

Bovine

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

LGB

Uniprot

P02754

Synonyms:

Beta-lactoglobulin (Beta-LG) (allergen Bos d 5), LGB


Immunogen:

Bovine milk Beta-lactoglobulin (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Bovine milk Beta-lactoglobulin antibody at 2 μ g/ml

Lane 1: Bovine milk Beta-lactoglobulin at 0.1mg/ml

Lane 2: Bovine milk Beta-lactoglobulin at 10 μ g/ml

Secondary

Goat polyclonal to rabbit IgG at 1/15000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa

.