

BTN3A1 Antibody, HRP conjugated



PACO25289

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in T-cell activation and in the adaptive immune response. Regulates the proliferation of activated T-cells. Regulates the release of cytokines and IFNG by activated T-cells. Mediates the response of T-cells toward infected and transformed cells that are characterized by high levels of phosphorylated metabolites, such as isopentenyl pyrophosphate.

Gene ID:

BTN3A1

Uniprot

O00481

Synonyms:

Butyrophilin subfamily 3 member A1 (CD antigen CD277), BTN3A1, BTF5

Immunogen:

Recombinant Human Butyrophilin subfamily 3 member A1 protein (30-254AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A