

# BTN3A1 Antibody, FITC conjugated



PACO25290

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Plays a role in T-cell activation and in the adaptive immune response. Regulates the proliferation of activated T-cells. Regulates the release of cytokines and IFNG by activated T-cells. Mediates the response of T-cells toward infected and transformed cells that are characterized by high levels of phosphorylated metabolites, such as isopentenyl pyrophosphate.

**Gene ID:**

BTN3A1

**Uniprot**

O00481

**Synonyms:**

Butyrophilin subfamily 3 member A1 (CD antigen CD277), BTN3A1, BTF5

**Immunogen:**

Recombinant Human Butyrophilin subfamily 3 member A1 protein (30-254AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A