

C19orf47 Antibody, Biotin conjugated



PACO25339

Product Information

Size:	Protein Background:
50ug	nucleoplasm.
Reactivity:	Gene ID:
Human	C19orf47
Source:	Uniprot
Rabbit	Q8N9M1
Isotype:	Synonyms:
IgG	Uncharacterized protein C19orf47, C19orf47
Applications:	Immunogen:
ELISA	Recombinant Human Uncharacterized protein C19orf47 protein (1-281AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A