

C1orf105 Antibody




PACO25340


Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: C1orf105
Source: Rabbit	Uniprot O95561
Isotype: IgG	Synonyms: Uncharacterized protein C1orf105, C1orf105
Applications: ELISA, IHC	Immunogen: Recombinant Human Uncharacterized protein C1orf105 protein (1-183AA).
Recommended dilutions: ELISA:1:2000-1:10000, IHC:1:20-1:200	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO25340 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO25340 at dilution of 1:100.