

Ovomucoid Antibody



PACO25388

Product Information

Size:

50ug

Reactivity:

Chicken

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

P01005

Synonyms:

Ovomucoid (Allergen Gal d I) (allergen Gal d 1)

Immunogen:

Chicken Ovomucoid (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A