

COX5A Antibody, FITC conjugated



PACO25402

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is the heme A-containing chain of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

Gene ID:

COX5A

Uniprot

P20674

Synonyms:

Cytochrome c oxidase subunit 5A, mitochondrial (Cytochrome c oxidase polypeptide Va), COX5A

Immunogen:

Recombinant Human Cytochrome c oxidase subunit 5A, mitochondrial protein (42-150AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A