

# COX5A Antibody, Biotin conjugated



PACO25403

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

This is the heme A-containing chain of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

**Gene ID:**

COX5A

**Uniprot**

P20674

**Synonyms:**

Cytochrome c oxidase subunit 5A, mitochondrial (Cytochrome c oxidase polypeptide Va), COX5A

**Immunogen:**

Recombinant Human Cytochrome c oxidase subunit 5A, mitochondrial protein (42-150AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A