

CA1 Antibody, Biotin conjugated



PACO25430

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Reversible hydration of carbon dioxide. Can hydrates cyanamide to urea.

Gene ID:

CA1

Uniprot

P00915

Synonyms:

Carbonic anhydrase 1 (EC 4.2.1.1) (Carbonate dehydratase I) (Carbonic anhydrase B) (CAB) (Carbonic anhydrase I) (CA-I), CA1

Immunogen:

Recombinant Human Carbonic anhydrase 1 protein (2-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A