

Casp12 Antibody



PACO25500

Product Information

Size:

50ug

Reactivity:

Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Involved in the activation cascade of caspases responsible for apoptosis execution.

Gene ID:

Casp12

Uniprot

O08736

Synonyms:

Caspase-12 (CASP-12) (EC 3.4.22. -), Casp12


Immunogen:

Recombinant Mouse Caspase-12 protein (1-92AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Casp12 antibody at 6 μ g/ml + LO2 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48 kDa

Observed band size: 48, 75 kDa