

Cd46 Antibody




PACO25600

Product Information

Size: 50ug	Protein Background:
Reactivity: Mouse	Gene ID: Cd46
Source: Rabbit	Uniprot O88174
Isotype: IgG	Synonyms: Membrane cofactor protein (CD antigen CD46), Cd46, Mcp
Applications: ELISA, WB, IHC	Immunogen: Recombinant Mouse Membrane cofactor protein (45-329AA).
Recommended dilutions: ELISA:1:2000-1:10000, WB:1:1000-1:5000, IHC:1:20-1:200	Storage: Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western Blot

Positive WB detected in: Mouse skeletal muscle tissue

All lanes: Cd46 antibody at 4 μ g/ml


Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 41, 34 kDa

Observed band size: 41 kDa

.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO25600 at dilution of 1:100.