

Cd46 Antibody, FITC conjugated



PACO25602

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Cd46

Uniprot

O88174

Synonyms:

Membrane cofactor protein (CD antigen CD46), Cd46, Mcp

Immunogen:

Recombinant Mouse Membrane cofactor protein (45-329AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A