

Cd63 Antibody, HRP conjugated



PACO25613

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play some role in signal transduction pathways.

Gene ID:

Cd63

Uniprot

P41731

Synonyms:

CD63 antigen (CD antigen CD63), Cd63

Immunogen:

Recombinant Mouse CD63 antigen protein (103-203AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A