

# Cd63 Antibody, Biotin conjugated



PACO25615

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play some role in signal transduction pathways.

**Gene ID:**

Cd63

**Uniprot**

P41731

**Synonyms:**

CD63 antigen (CD antigen CD63), Cd63

**Immunogen:**

Recombinant Mouse CD63 antigen protein (103-203AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A