

# Cd74 Antibody, Biotin conjugated



PACO25630

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Plays a critical role in MHC class II antigen processing by stabilizing peptide-free class II alpha/beta heterodimers in a complex soon after their synthesis and directing transport of the complex from the endoplasmic reticulum to compartments where peptide loading of class II takes place.

**Gene ID:**

Cd74

**Uniprot**

P04441

**Synonyms:**

H-2 class II histocompatibility antigen gamma chain (Ia antigen-associated invariant chain) (Ii) (MHC class II-associated invariant chain) (CD antigen CD74), Cd74, Ii

**Immunogen:**

Recombinant Mouse H-2 class II histocompatibility antigen &gamma; chain protein (56-279AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A