

Cd82 Antibody



PACO25636

Product Information

Size:

50ug

Reactivity:

Mouse, Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Gene ID:

Cd82

Uniprot

P40237

Synonyms:

CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1 homolog) (CD antigen CD82), Cd82, Kai1


Immunogen:

Recombinant Mouse CD82 antigen protein (111-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Cd82 antibody at 2 μ g/ml

Lane 1: Jurkat whole cell lysate

Lane 2: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 30 kDa

Observed band size: 30 kDa

.