

CDC42EP2 Antibody, HRP conjugated



PACO25665

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probably involved in the organization of the actin cytoskeleton. May act downstream of CDC42 to induce actin filament assembly leading to cell shape changes. Induces pseudopodia formation in fibroblasts in a CDC42-dependent manner.

Gene ID:

CDC42EP2

Uniprot

O14613

Synonyms:

Cdc42 effector protein 2 (Binder of Rho GTPases 1), CDC42EP2, BORG1 CEP2

Immunogen:

Recombinant Human Cdc42 effector protein 2 protein (9-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A