

CDK2AP1 Antibody



PACO25672

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

specific inhibitor of the cell-cycle kinase CDK2.

Gene ID:

CDK2AP1

Uniprot

O14519

Synonyms:

Cyclin-dependent kinase 2-associated protein 1 (CDK2-associated protein 1) (Deleted in oral cancer 1) (DOC-1) (Putative oral cancer suppressor), CDK2AP1, CDKAP1 DOC1


Immunogen:

Recombinant Human Cyclin-dependent kinase 2-associated protein 1 protein (1-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO25672 at dilution of 1:100.