

CDK2AP1 Antibody, FITC conjugated



PACO25674

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

specific inhibitor of the cell-cycle kinase CDK2.

Gene ID:

CDK2AP1

Uniprot

O14519

Synonyms:

Cyclin-dependent kinase 2-associated protein 1 (CDK2-associated protein 1) (Deleted in oral cancer 1) (DOC-1) (Putative oral cancer suppressor), CDK2AP1, CDKAP1 DOC1

Immunogen:

Recombinant Human Cyclin-dependent kinase 2-associated protein 1 protein (1-115AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A