

# CDK2AP1 Antibody, Biotin conjugated



PACO25675

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

specific inhibitor of the cell-cycle kinase CDK2.

**Gene ID:**

CDK2AP1

**Uniprot**

O14519

**Synonyms:**

Cyclin-dependent kinase 2-associated protein 1 (CDK2-associated protein 1) (Deleted in oral cancer 1) (DOC-1) (Putative oral cancer suppressor), CDK2AP1, CDKAP1 DOC1

**Immunogen:**

Recombinant Human Cyclin-dependent kinase 2-associated protein 1 protein (1-115AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A