

CDK2AP2 Antibody



PACO25676

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

CDK2AP2

Uniprot

O75956

Synonyms:

Cyclin-dependent kinase 2-associated protein 2 (CDK2-associated protein 2) (DOC-1-related protein) (DOC-1R), CDK2AP2, DOC1R


Immunogen:

Recombinant Human Cyclin-dependent kinase 2-associated protein 2 protein (1-126AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO25676 at dilution of 1:100.