

CHAD Antibody, Biotin conjugated



PACO25732

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Promotes attachment of chondrocytes, fibroblasts, and osteoblasts. This binding is mediated (at least for chondrocytes and fibroblasts) by the integrin alpha(2)beta(1). May play an important role in the regulation of chondrocyte growth and proliferation.

Gene ID:

CHAD

Uniprot

O15335

Synonyms:

Chondroadherin (Cartilage leucine-rich protein), CHAD, SLRR4A

Immunogen:

Recombinant Human Chondroadherin protein (23-359AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A