

ODC1 Antibody, HRP conjugated



PACO25746

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

L-ornithine = putrescine + CO₂. Inhibited by S-nitrosylation.

Gene ID:

ODC1

Uniprot

P11926

Synonyms:

Ornithine decarboxylase (ODC) (EC 4.1.1.17), ODC1

Immunogen:

Recombinant Human Ornithine decarboxylase protein (1-461AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A