

# ODC1 Antibody, FITC conjugated



PACO25747

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

L-ornithine = putrescine + CO<sub>2</sub>. Inhibited by S-nitrosylation.

**Gene ID:**

ODC1

**Uniprot**

P11926

**Synonyms:**

Ornithine decarboxylase (ODC) (EC 4.1.1.17), ODC1

**Immunogen:**

Recombinant Human Ornithine decarboxylase protein (1-461AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A