

# CHRM3 Antibody, FITC conjugated



PACO25755

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

The muscarinic acetylcholine receptor mediates various cellular responses, including inhibition of adenylate cyclase, breakdown of phosphoinositides and modulation of potassium channels through the action of G proteins. Primary transducing effect is Pi turnover.

**Gene ID:**

CHRM3

**Uniprot**

P20309

**Synonyms:**

Muscarinic acetylcholine receptor M3, CHRM3

**Immunogen:**

Recombinant Human Muscarinic acetylcholine receptor M3 protein (253-492AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A