

# CKS2 Antibody, HRP conjugated



PACO25770

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.

**Gene ID:**

CKS2

**Uniprot**

P33552

**Synonyms:**

Cyclin-dependent kinases regulatory subunit 2 (CKS-2), CKS2

**Immunogen:**

Recombinant Human Cyclin-dependent kinases regulatory subunit 2 protein (2-79AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A