

CKS2 Antibody, Biotin conjugated



PACO25772

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function.

Gene ID:

CKS2

Uniprot

P33552

Synonyms:

Cyclin-dependent kinases regulatory subunit 2 (CKS-2), CKS2

Immunogen:

Recombinant Human Cyclin-dependent kinases regulatory subunit 2 protein (2-79AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A